Porous Power

Poröse Power

By developing CoralPor™ 1000 and CoralPor™ 2000, researchers at SCHOTT have introduced a porous glass product that has applications in a variety of sectors.

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Photo: foto: schott/H.-J. Schulz

THOMAS H. LOEWE

They are inconspicuous but very important: in labs around the world, reference electrodes detect single particles among trillions. As such, they are inevitably part of almost all electrochemical measurements. Tiny porous plugs are often used to separate the reference fluid inside the reference electrodes from the sample. For this, SCHOTT now offers CoralPor™ 1000, a nanoporous glass solution that improves measurements and is even suitable for applications far beyond the laboratory.

Seen under a scanning electron microscope, CoralPor™ resembles the tender structure of name-giving coral skeletons. At high magnification, countless pores and channels are revealed. The interconnected coral-like structure is created during a special production process, in which borosilicate glass is melted, rapidly quenched and then phase separated. This process triggers separation of the individual components in the glass, resulting in a sodium borate phase and a silicate phase. At the end of the process, special chemicals are used to dissolve the sodium borate phase, which yields an interconnected microstructure with open porosity. The innovative glass is currently available in two product variants: CoralPor™ 1000 is typically employed in applications such as reference electrode junctions, desiccants, coatings, or as components in devices for the medical, aerospace, defense, and oil and gas industries. CoralPor™ 2000 exhibits larger pore sizes (40–300 nm) and is best suited for use in separations or as a synthesis substrate. Both CoralPor™ products were developed

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under the leadership of Dr. William James, Supervisor of Materials Development at SCHOTT North America. "Our goal was simple: to produce a porous glass that offered high chemical and mechanical stability," says James.

One example that clearly demonstrates CoralPor™ 1000's extreme stability: "CoralPor™ can even be utilized as a key element for heat-resistant coatings on spacecraft," says Ed Hart, Senior Manager of Market Development and Innovation at SCHOTT North America. Since the material consists of 95-97 percent silica, it also has a low thermal expansion coefficient. This renders the glass extremely robust for thermal protection applications and ensures that it expands minimally – a major advantage at high temperatures, which makes CoralPor™ 1000 a "universally" applicable material in the truest sense of the word. But Hart also emphasizes another main advantage: "CoralPor™ 1000 has lots of tiny pores on the order of 4 nanometers. This makes the product ideal for reference electrodes and other lab applications." In fact, lab applications were the actual inspiration for developing CoralPor™ 1000.

For many years nanoporous glass was manufactured by one supplier. But the recent discontinuation of the product caused major concern in the scientific community. Between 2012 and 2013 alone, more than 70 publications in various scientific journals explicitly reported using reference electrodes with glass plugs sold under the now discontinued brand. "We were quite worried," says Philippe Buhlmann, Professor of Chemistry and Chemical Physics at the University of Minnesota. "But now, many members of the community are relieved to see that SCHOTT is stepping in with its CoralPor™ product line," Buhlmann adds. CoralPor™ has also gained trust among users of chromatography applications – a

The nanoporous glass CoralPor™ (photo p. 36) has a broad spectrum of applications, e.g. in the area of electrochemistry, chromatography (see below) and for heat-resistant coatings on spacecrafts (right).

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
Photo Foto : Thinkstock

chemical separation method used in the purification of biotechnology products. For this process, researchers require the separation media to withstand up to pH 14 to survive the cleaning process between chromatography cycles. "Up to now, conventional chromatography glass-based media could only remain chemically

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Dr. Jonathon O. Howell, a Senior Scientist at Emergent Instruments, is an expert in electrochemistry – and provided the inspiration for the development of CoralPor™ 1000.

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solutions: How exactly do reference electrodes work?

Howell: Reference electrodes serve as a point of reference when measuring potential to define the pH value of a chemical solution for example. A reference electrode is a unit separated by glass that is filled with a standard salt solution. Chemical reactions take place within it and produce electric potential. Dipping the reference electrode in a measuring solution elicits a potential difference. This can then be used to calculate the measurement value, for example the pH value.

solutions: Why do you need porous glass for that?

Howell: Potential difference can only exist if the solution on the inside of the glass is in direct contact with the measuring solution on the outside. Porous glass frits, for example, allow a minimal amount of ions to flow through their pores and channels. However, at the same time these openings in the glass are so small that neither the measuring nor the standard solutions are contaminated. The researchers at SCHOTT involved me very closely in the development of CoralPor™ and it has proven to be an ideal material for our glass frits.

solutions: In your opinion, what are the advantages of working with SCHOTT?

Howell: CoralPor™ is always available quickly and in sufficient quantities. SCHOTT is also very open to customer feedback and excels in adapting to the wishes of its customers. After all, that was how CoralPor™ 1000 came about – as the result of an e-mail I sent them.

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Photo Foto : Emergent Instruments

stable in a range of pH 2 to pH 8,” says James. The specific surface area of CoralPor™ delivers an advantage by maximizing the amount of target molecules separated per cycle.

“The great thing about CoralPor™ is that we can adapt so much to the customer’s specific demands,” says Hart. SCHOTT researchers can tailor CoralPor™ to a wide range of geometries, surface areas and pore sizes. This flexibility also shows that CoralPor™ has a broad spectrum of possible applications remaining to be explored: “One direction we’re interested in is to further increase the surface area of CoralPor™ in turn increasing the material’s capability as an absorbent material and synthesis substrate,” says James. Heat-resistant coatings for turbines and engines are another possibility. All this shows: the experts at SCHOTT won’t be running out of ideas any time soon.

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
1.

MELT:

Silica-based sand and additional components are mixed and melted.

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
2.

PHASE SEPARATION:

The quenched glass from the melt is put into an annealing lehr imparting the interconnected microstructure through phase separation.

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3.

CHEMICAL LEACHING:

The phase separated glass is put into a bath to chemically leach the soluble sodium borate phase and achieve open porosity.

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4.

FORMATTING:

After leaching, different shapes are fabricated to meet customer specifications.

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