




Footswitch



Converter Cable



Polarizing filters



Diffusor and protection window for ringlight

Accessories for VisiLED

For best performance of LED illumination

Footswitch

- Enables the user to remotely switch between stored memory positions 1 through 4 to quickly change the stored illumination conditions.
- Allows the user to quickly and easily change to various illumination types by allowing total hand control under the microscope.

RS-232 USB 1.1 Converter Cable

- Connects the serial plug of the MC 1500 controller to the USB plug of a computer.
- Full USB 1.1 functionality will be maintained.
- Driver is included.

Electronic Accessories

Description	Part No.	Length
Footswitch	400 020	
RS-232 USB 1.1 Converter Cable	400 030	1.5 meter (59")

Polarization filters

- Polarization system for brightfield ringlights: turnable polarizer with adapting thread and sliding filter as analyzer.
- Turnable polarization filter attachment for transmitted brightfield light stage.

Polarization Filters

	Part No.	Mounting
Polarization filter system for ringlights	400 550	Polarizer M110 x 1 Analyzer of 66 mm sliding drawer
Polarization filter attachment for transmitted light stage	158 500	Ø 84 mm
Analyzer	158 505	M49 x 0.75
Analyzer	158 510	M52 x 1

Accessories for ringlights

- Diffusor with adapting thread (M110 x 1) to avoid reflections of LEDs.
- Borofloat® glass protection window with adapting thread (M110 x1). Protects LEDs and microscope objective in oily and dusty environments.

Accessories for ringlights

	Part No.	Mounting
Diffusor	400 570	M110 x 1
Protection window	400 560	M110 x 1

*Product improvements may result in
specification or feature changes without notice.*

Lighting and Imaging
SCHOTT North America, Inc.
122 Charlton Street
Southbridge, MA 01550-1960
USA

Phone: +1 508/765-9744
Fax: +1 508/764-6273
lightingimaging@us.schott.com
www.us.schott.com/lightingimaging

For your local sales contact please visit
www.us.schott.com/lightingimaging

SCHOTT
glass made of ideas