

Electronic Accessories for KL 2500 LCD

The KL 2500 LCD has a built in 9-pin analog plug, which enables the connection of electronic accessories.



Remote control, P/N 258 402



Cable extension, P/N 258 401



Y-piece, P/N 258 400



Tilt switch, P/N 258 403

Features

- Remote control: analog control pad for remote setting of intensity. Switches the light sources also in “on” or “standby”-modes.
- Cable extension: 10 meter extension for remote control or tilt switch.
- Y-piece: allows simultaneous use of remote control and tilt switch.
- Tilt switch: enables the light source to be mounted onto a swivel-arm, i. e. for use with surgical or diagnostic microscopes. It switches from nonworking / standby to working position: i.e. when the arm is moved upwards the tilt switch switches the light source off. When the arm is moved down it switches it on again. The standby mode is displayed on the LCD display of the light source.

Electronic accessories for KL 2500 LCD	Part. No.	Dimension
Remote control	258 402	3 m
Cable extension	258 401	10 m
Y-piece	258 400	
Tilt switch	258 403	

All specifications are subject to change without prior notice. This datasheet or any extracts thereof may only be used in other publications with express permission of SCHOTT. © SCHOTT North America, Inc.

SCHOTT North America, Inc.
Lighting and Imaging
 6862 Santa Teresa Blvd.
 San Jose, CA 95119 USA
 phone: +1-408-363-2100
 fax: +1-408-363-9980
 microscopy@us.schott.com
 www.us.schott.com/lightingimaging

