

# Accessories for VisiLED

For best performance of LED illumination




Footswitch



Converter Cable



Polarizing filters



Diffusor and protection window for ringlight

## Footswitch

- Enables the user to remotely switch between stored memory positions M1 ... M4 to quickly change the stored illumination conditions.

## RS-232 USB 1.1 Converter Cable

- Connects the serial plug of the MC 1500 controller to the USB plug of a computer.
- Full USB 1.1 functionality will be maintained.
- Driver is included.

## Polarization filters

- Polarization system for brightfield ringlights: turnable polarizer with adapting thread and sliding filter as analyzer.
- Turnable polarization filter attachment for transmitted brightfield light stage.

## Accessories for ringlights

- Diffusor with adapting thread to avoid reflections of LEDs.
- Borofloat® glass protection window with adapting thread, to protect LEDs and microscope objective in oily and dusty environments.

Polarization Filters		
Description	Part. No.	Mounting
Polarization filter set for S80-55 brightfield ringlight.	400 550	Polarizer M110 x 1 Analyzer for 66 mm lens
Polarization filter set for S80-65N brightfield ringlight.	400 665	Polarizer M130 x 1 Analyzer for 80 mm lens
Polarization filter set for S40-75 brightfield ringlight.	400 575	Polarizer M90 x 1 Analyzer for 66 mm lens
Polarization filter attachment for transmitted light stage	158 500	Ø 84 mm
Analyzer	158 505	M49 x 0.75
Analyzer	158 510	M52 x 1

Accessories for ringlights			Electronic accessories		
	Part. No.	Mounting		Part. No.	Length
Diffusor (not for S80-65N, S40-75)	400 570	M110 x 1	Footswitch, momentary	400 020	
Diffusor for S40-75 brightfield ringlight	400 580	M90 x 1		RS-232 USB 1.1 Converter Cable	400 030
Protection window	400 560	M110 x 1			

All specifications are subject to change without prior notice. This datasheet or any extracts thereof may only be used in other publications with express permission of SCHOTT. © SCHOTT North America, Inc.

SCHOTT North America, Inc.  
**Lighting and Imaging**  
 122 Charlton Street  
 Southbridge, MA 01550 USA  
 Phone: +1-508-765-3211  
 Fax: +1-508-764-6273  
 lightingimaging@us.schott.com  
 www.us.schott.com/lightingimaging



Version 03.2018