



Ringlight to Microscope Adapters


A08606




A08607




A08608




A08609



A08626



A08628



Ringlight Adapters

Note:

Adapters are machined to order. Pricing and availability subject to market conditions.

Drawings are available.

Part No.	Ringlight	Microscope	Mounting
A08606	A08600	B&L Zoom 4	Screws into Objective
A08607	A08625	Leica SZ4/SZ6 with auxiliary lens	Three Thumbscrews
A08608	A08625	Leica SZ4/SZ6 without auxiliary lens	Screws into objective
A08609	A08630	Leica SZ7	#6 - 32 Hex Set Screw
A08626	A08625	Olympus SZ	Screw-on Retaining Ring
A08628	A08600, A08620	SMZ 645/660	Screws on Objective

Dimensions in () are in mm

All specifications are subject to change without prior notice. This datasheet or any extracts thereof may only be used in other publications with express permission of SCHOTT. © SCHOTT North America, Inc.

Lighting and Imaging
SCHOTT North America, Inc.
 122 Charlton Street
 Southbridge, MA 01550
 USA

Phone: +1 (508) 765-9744
 Fax: +1 (508) 765-1299
 lightingimaging@us.schott.com
 www.us.schott.com/lightingimaging

