

Controllers and Power Supply for VisiLED Series



MC750
Basic functions



MC1000
Compact and ergonomic



MC1500
Multifunction center



CEE 7/7 plug



NEMA 5-15 plug



BS 1363 plug



SEV1011 plug

	MC750	MC1000	MC1500
Dimensions	73.5 mm x 110 mm x 40 mm (2.9" x 4.3" x 1.58")	90 mm x 152 mm x 25 mm (3.54" x 5.98" x 0.98")	153 mm x 110 mm x 40 mm (6" x 4.3" x 1.58")
Connectors on backside	<ul style="list-style-type: none"> • 21 - 24V DC Connector • 1 Lighthouse Connector 	<ul style="list-style-type: none"> • 21V DC Connector • 1 Lighthouse Connector • ESD plug 	<ul style="list-style-type: none"> • 21 - 24V DC Connector • RS232 Interface • 2 Lighthouse Connectors • 1 Footswitch • 1 Flash
Features	<ul style="list-style-type: none"> • Continuous dimming via potentiometer 	<ul style="list-style-type: none"> • Easy to use. • 4 present segment modes. • Rotation to both directions. • Continuous dimming via potentiometer. 	<ul style="list-style-type: none"> • 5 preset segment modes plus rotation for optimum contrasting. • Smooth dimming with overdrive position. • 4 memory positions to save custom lighting modes. • External Trigger: to start single strobe. • Flash mode: to overdrive continuous set intensity with a short flash for taking photographs. • Strobe mode: to utilize internal strobe function with frequencies between 15 Hz – 5 KHz. • RS 232 Computer interface (USB 1.1 optional) to set parameters via computer. • Foot switch between memory positions (optional).




Power Supplies	
Input 100–240V AC, 50/60 Hz; Output 21V DC/1,92A, UL certification and CE marking	400 050
Power Cord IEC-C13 to CEE 7/7 (European Plug), 1.8 m	400 051
Power Cord IEC-C13 to NEMA 5-15 (US Plug), 1.8 m	400 052
Power Cord IEC-C13 to BS 1363 (UK Plug), 1.8 m	400 053
Power Cord IEC-C13 to SEV1011 (Swiss Plug), 1.8 m	400 054

Controller Model	Part. No.
MC750	400 010
MC1000	400 040
MC1500	400 000




Version 11.2018






Controllers and Power Supply for VisiLEDs MC Series



VisiLED controller MC750



VisiLED controller MC1000



VisiLED controller MC1500

All specifications are subject to change without prior notice. This datasheet or any extracts thereof may only be used in other publications with express permission of SCHOTT. © SCHOTT North America, Inc.

Lighting and Imaging
SCHOTT North America, Inc.
 122 Charlton Street
 Southbridge, MA 01550 USA
 Phone: +1-508-765-3211
 Fax: +1-508-764-6273
 lightingimaging@us.schott.com
 www.us.schott.com/lightingimaging

